

# **Faculty Research Lab**

## **(Physics & Chemistry)**



Section	Equipment	Page
1	Bruker D6 Phaser X-Ray Diffractometer (XRD)	3
2	Agilent Gas Chromatography–Mass Spectrometry (GC-MS)	4
3	Renishaw inVia Raman Microscope	5
4	Tabletop Atomic Force Microscope (AFM)	6
5	KinetAsyst Stopped-Flow Spectrometer	7
6	Thermo Scientific Nicolet 380 ATR-FTIR (Teaching FTIR)	8
7	Thermo Scientific Nicolet 4700 FT-IR Spectrophotometer (Research FTIR)	9
8	Agilent Cary 60 UV-Vis Spectrophotometer	10
9	Agilent Cary 8454 UV-Vis Spectrophotometer	11
10	Anasazi Eft-60 Nuclear Magnetic Resonance (NMR) Spectrometer	12
11	Buehler MetaServ250 Metallographic Grinder/Polisher	13
12	Film Sense Multi-Wavelength Ellipsometer	14
13	Glancing Angle Deposition (GLAD) System	15
14	Harrick PDC Plasma Cleaner	16
15	ISCPS Corp nGauge Atomic Force Microscope	17
16	Perkin Elmer Analyst 100 Atomic Absorption Spectrophotometer	18
17	StellarNet Fluor-System Fluorescence Spectrometer	19
18	Thermo Scientific Barnstead Smart2Pure Water Purifier	20
19	Thermo Scientific LTQ XL Mass Spectrometer	21
20	Varian CP-3380 Gas Chromatograph	22
21	Mbraun Labstar Glovebox	23
22	Dionex UltiMate 3000 UHPLC System	24
23	Pine Bipotentiostat (Model AFCBP1)	25
24	Muffle Furnace KSL-1200X	26
25	Heidolph Rotary Evaporator	27
26	Nima Technology Langmuir-Blodgett Trough	28
27	Horiba DeltaPro Lifetime Fluorimeter	29

# 1. Bruker D6 Phaser X-Ray Diffractometer (XRD)



## Description

The Bruker D6 Phaser is a compact benchtop X-ray diffractometer used for crystallographic analysis and phase identification of solid materials. The system enables researchers to determine crystal structure, phase composition, and lattice parameters for a wide range of inorganic and organic materials.

## Capabilities

- Copper X-ray source with high stability
- Silicon strip detector for high-resolution diffraction data
- Angular measurement range from  $-3^{\circ}$  to  $152^{\circ}$  ( $2\theta$ )
- Automated sample alignment and measurement
- Integrated analysis software for diffraction pattern interpretation

## Applications

- Crystal structure identification
- Phase analysis of engineering and chemical materials
- Thin film and powder diffraction studies
- Materials science and solid-state chemistry research

## 2. Agilent Gas Chromatography–Mass Spectrometry (GC-MS)



### Description

The Agilent GC-MS system combines gas chromatography with mass spectrometry to separate, identify, and quantify volatile and semi-volatile chemical compounds in complex mixtures.

### Capabilities

- High-resolution gas chromatographic separation
- Mass spectral detection for compound identification
- Automated sample injection
- Integrated spectral libraries for compound analysis

### Applications

- Environmental contaminant analysis
- Pharmaceutical and biochemical analysis
- Chemical composition identification
- Trace compound detection in research samples

### 3. Renishaw inVia Raman Microscope



#### Description

The Renishaw inVia Raman Microscope is a high-resolution Raman spectroscopy system used for molecular and structural characterization of materials at microscopic scales.

#### Capabilities

- High-sensitivity Raman spectral detection
- Confocal microscope integration for spatial analysis
- Multiple excitation laser options
- High spectral resolution and rapid data acquisition

#### Applications

- Molecular structure identification
- Chemical composition analysis
- Semiconductor and nanomaterial characterization
- Pharmaceutical and polymer research

## 4. Tabletop Atomic Force Microscope (AFM)



### Description

The Tabletop Atomic Force Microscope is a scanning probe microscope used to measure surface topography and nanoscale features of materials with extremely high resolution.

### Capabilities

- Nanometer-scale surface imaging
- Contact and non-contact scanning modes
- High-precision cantilever sensors
- Digital surface mapping and analysis

### Applications

- Nanomaterials characterization
- Surface morphology analysis
- Thin-film research
- Semiconductor and materials science studies

## 5. KinetAsyst Stopped-Flow Spectrometer



### Description

The KinetAsyst Stopped-Flow Spectrometer is designed to measure extremely fast chemical reactions by rapidly mixing reactants and monitoring changes in absorbance or fluorescence over very short time intervals.

### Capabilities

- Millisecond reaction monitoring
- Rapid mixing of chemical reactants
- Dual optical detection channels
- Automated kinetic data acquisition

### Applications

- Enzyme kinetics research
- Biochemical reaction analysis
- Pharmaceutical reaction studies
- Chemical kinetics experimentation

## 6. Thermo Scientific Nicolet 380 ATR-FTIR (Teaching FTIR)



### Description

The Thermo Scientific Nicolet 380 ATR-FTIR Spectrometer is an infrared spectroscopy instrument used primarily for instructional laboratories to analyze molecular vibrations and identify chemical functional groups.

### Capabilities

- ATR sampling accessory for solid and liquid samples
- Mid-infrared spectral analysis
- Rapid spectral acquisition
- User-friendly teaching interface

### Applications

- Organic compound identification
- Functional group analysis
- Chemistry laboratory instruction
- Polymer and material characterization

## 7. Thermo Scientific Nicolet 4700 FT-IR Spectrophotometer (Research FTIR)



### Description

The Thermo Scientific Nicolet 4700 FT-IR Spectrophotometer is a high-performance infrared spectroscopy system used for advanced molecular and structural analysis in research laboratories.

### Capabilities

- High-resolution FT-IR spectral analysis
- Broad infrared wavelength detection range
- Modular accessory compatibility
- Advanced spectral processing software

### Applications

- Molecular structure analysis
- Chemical reaction monitoring
- Polymer characterization
- Advanced materials research

## 8. Agilent Cary 60 UV-Vis Spectrophotometer



### Description

The Agilent Cary 60 UV-Vis Spectrophotometer measures absorbance and transmittance of liquid samples across ultraviolet and visible wavelengths for quantitative chemical analysis.

### Capabilities

- Wavelength range from 190 to 1100 nm
- Xenon flash lamp light source
- Rapid data acquisition
- High sensitivity for trace analysis

### Applications

- DNA and protein quantification
- Chemical kinetics experiments
- Environmental analysis
- Pharmaceutical research

## 9. Agilent Cary 8454 UV-Vis Spectrophotometer



### Description

The Agilent Cary 8454 UV-Vis Spectrophotometer is a diode-array spectrophotometer designed for rapid spectral analysis and monitoring of chemical reactions.

### Capabilities

- Diode-array detection system
- Simultaneous multi-wavelength measurement
- Rapid full-spectrum acquisition
- Automated data analysis software

### Applications

- Reaction monitoring
- Quantitative chemical analysis
- Spectroscopic research studies
- Environmental and pharmaceutical testing

# 10. Anasazi Eft-60 Nuclear Magnetic Resonance (NMR) Spectrometer



## Description

The Anasazi Eft-60 NMR Spectrometer is an analytical instrument used to determine the molecular structure and chemical composition of organic compounds through nuclear magnetic resonance spectroscopy.

## Capabilities

- High-field NMR analysis
- Digital spectral acquisition and processing
- Automated sample analysis
- High-precision frequency stability

## Applications

- Molecular structure determination
- Organic compound identification
- Chemical synthesis verification
- Pharmaceutical and biochemical research

# 11. Buehler MetaServ250 Metallographic Grinder/Polisher



## Description

The Buehler MetaServ250 is a metallographic preparation system used to grind and polish material samples to produce smooth surfaces suitable for microscopic and structural analysis.

## Capabilities

- Variable speed grinding and polishing
- Automated polishing controls
- Multi-sample preparation capability
- Durable industrial design

## Applications

- Metallurgical sample preparation
- Microstructure analysis
- Materials science research
- Quality control testing

## 12. Film Sense Multi-Wavelength Ellipsometer



### Description

The Film Sense Multi-Wavelength Ellipsometer is an optical instrument used to measure thin-film thickness and optical properties of materials using polarized light analysis.

### Capabilities

- Multi-wavelength ellipsometry measurements
- High-precision thin film thickness measurement
- Automated optical property analysis
- Computer-controlled measurement software

### Applications

- Thin film characterization
- Semiconductor research
- Nanotechnology studies
- Materials science analysis

# 13. Glancing Angle Deposition (GLAD) System



## Description

The Glancing Angle Deposition System is a thin-film fabrication system used to create nanostructured coatings through controlled vapor deposition at oblique angles.

## Capabilities

- Controlled deposition angles
- Thin-film growth monitoring
- Vacuum deposition environment
- Adjustable substrate positioning

## Applications

- Nanostructured material fabrication
- Thin-film optical coatings
- Semiconductor device research
- Advanced materials development

# 14. Harrick PDC Plasma Cleaner



## Description

The Harrick PDC Plasma Cleaner is a plasma surface treatment system used to clean and activate surfaces prior to thin-film deposition or bonding processes.


## Capabilities

- Radio-frequency plasma generation
- Surface contamination removal
- Adjustable plasma exposure parameters
- Compact benchtop design

## Applications

- Surface cleaning for microscopy samples
- Semiconductor device preparation
- Surface modification research
- Materials science experiments

# 15. ISCPs Corp nGauge Atomic Force Microscope



## Description

The ISCPs nGauge Atomic Force Microscope is a high-precision scanning probe microscope used to characterize nanoscale surface structures and material properties.

## Capabilities

- Nanometer-scale imaging resolution
- Multiple scanning modes
- High-precision cantilever detection
- Digital surface reconstruction software

## Applications

- Nanotechnology research
- Surface roughness analysis
- Thin film characterization
- Semiconductor device analysis

# 16. Varian CP-3380 Gas Chromatograph



## Description

The Varian CP-3380 Gas Chromatograph is an analytical instrument used to separate and analyze volatile compounds in complex mixtures. It employs gas chromatography to separate components based on their interactions with a stationary phase inside a capillary column.

## Capabilities

- Compatible with multiple detectors such as FID, TCD, and ECD
- High-resolution capillary column separation
- Automated temperature programming
- Computer-controlled data acquisition and analysis

## Applications

- Environmental contaminant analysis
- Pharmaceutical compound identification
- Petrochemical analysis
- Chemical mixture characterization

# 17. Mbraun Labstar Glovebox



## Description

The Mbraun Labstar Glovebox is a controlled-atmosphere enclosure designed for handling air-sensitive materials in an oxygen- and moisture-free environment.

## Capabilities

- Integrated gas purification system
- Low oxygen and moisture atmosphere control
- Dual glove ports for sample manipulation
- Real-time monitoring of environmental conditions

## Applications

- Handling air-sensitive chemicals
- Materials synthesis and preparation
- Battery and semiconductor research
- Controlled environment experiments

# 18. Dionex UltiMate 3000 UHPLC System



## Description

The Dionex UltiMate 3000 Ultra-High-Performance Liquid Chromatography (UHPLC) system is used for separating, identifying, and quantifying compounds in liquid samples with extremely high resolution.

## Capabilities

- High-pressure liquid chromatography operation up to ~1000 bar
- Automated sample injection and solvent delivery
- High-resolution separation columns
- Integrated software for data processing and analysis

## Applications

- Pharmaceutical compound analysis
- Environmental sample testing
- Biomolecule separation and analysis
- Chemical reaction monitoring

# 19. Pine Bipotentiostat (Model AFCBP1)



## Description

The Pine Bipotentiostat is an electrochemical instrument used to control and measure voltage and current in electrochemical experiments involving two working electrodes.


## Capabilities

- Potentiostatic and galvanostatic control modes
- High-precision voltage and current measurement
- Compatible with multiple electrochemical cells
- Integrated electrochemical data analysis software

## Applications

- Electrochemical reaction analysis
- Corrosion studies
- Battery and fuel cell research
- Electrochemical sensor development

## 20. Table Top Atomic Force Microscope



### Description

The Table Top Atomic Force Microscope is a compact scanning probe microscope used to analyze surface topography and nanoscale properties of materials.

### Capabilities

- Nanometer-scale surface imaging
- Contact and tapping scanning modes
- Digital surface mapping software
- Compact benchtop instrument design

### Applications

- Surface morphology characterization
- Nanomaterials research
- Thin film analysis
- Materials science and semiconductor research

# 21. HC-E1 Sumitomo Cryostat



## Description

The HC-E1 Sumitomo Cryostat is a cryogenic cooling system used to maintain extremely low temperatures for experiments requiring cryogenic environments.

## Capabilities

- Temperature control down to approximately 1.5 K
- Stable cryogenic cooling environment
- Compatible with various sample holders
- Automated temperature monitoring and control

## Applications

- Superconductivity research
- Low-temperature physics experiments
- Quantum materials research
- Cryogenic material characterization

## 22. Muffle Furnace – KSL-1200X



### Description

The KSL-1200X Muffle Furnace is a high-temperature laboratory furnace used for thermal processing, heat treatment, and material sintering.

### Capabilities

- Maximum temperature up to approximately 1200 °C
- Programmable temperature control
- Uniform heating chamber
- Durable high-temperature insulation

### Applications

- Heat treatment of materials
- Ceramic sintering
- Thermal decomposition experiments
- Materials science research

## 23. Heidolph Rotary Evaporator



### Description

The Heidolph Rotary Evaporator is a laboratory instrument used to remove solvents from samples through controlled evaporation under reduced pressure.


### Capabilities

- Adjustable rotation speed for efficient evaporation
- Vacuum-assisted solvent removal
- Controlled heating bath for temperature regulation
- Compatible with multiple glassware configurations

### Applications

- Solvent evaporation and concentration
- Chemical synthesis purification
- Natural product extraction
- Pharmaceutical sample preparation

## 24. Nima Technology Langmuir-Blodgett Trough



### Description

The Nima Technology Langmuir-Blodgett Trough is an instrument used to create and study thin molecular films at the air-water interface.

### Capabilities

- Controlled monolayer film deposition
- Precision barrier movement for surface pressure control
- High-accuracy film transfer onto substrates
- Integrated software for experimental control

### Applications

- Thin-film fabrication
- Nanomaterial research
- Surface chemistry studies
- Sensor and device development

# 25. Horiba DeltaPro Lifetime Fluorimeter



## Description

The Horiba DeltaPro Lifetime Fluorimeter is an advanced fluorescence spectroscopy system used to measure fluorescence lifetimes and molecular interactions.

## Capabilities

- Time-correlated single photon counting (TCSPC) detection
- Picosecond to nanosecond lifetime measurement
- High sensitivity fluorescence detection
- Modular optical configuration

## Applications

- Fluorescence lifetime analysis
- Molecular interaction studies
- Biophysics and nanomaterial research
- Photophysics and photochemistry experiments

## 26. Thermo Scientific Barnstead Smart2Pure Water Purifier



### Description

The Thermo Scientific Barnstead Smart2Pure Water Purifier provides high-purity laboratory water for analytical experiments and sample preparation.

### Capabilities

- Produces Type 1 and Type 2 purified water
- Integrated UV purification system
- Ion-exchange and filtration purification stages
- Digital monitoring of water quality

### Applications

- Preparation of analytical samples
- Laboratory reagent preparation
- Instrument feed water supply
- Chemical and biological experiments

## 27. StellarNet Fluor-System Fluorescence Spectrometer



### Description

The StellarNet Fluor-System Fluorescence Spectrometer is a compact spectroscopic instrument used to measure fluorescence emission from chemical and biological samples.

### Capabilities

- Broad excitation and emission wavelength detection
- High sensitivity fluorescence measurement
- Portable spectrometer design
- Integrated spectral analysis software

### Applications

- Fluorescence spectroscopy analysis
- Environmental monitoring
- Biochemical research
- Quantitative fluorescence measurements